





Dynamic PET Study + masks or TACs

LOADING DATA



MODEL SELECTION & PREPROCESSING

Kinetic model selection + input parameters for the preprocessing



Results
SRTM preprocessing:
R1: 1.2503
k2a: 0.070609
k2: 0.35891
BPnd: 4.083

Goodness of fit:
NMSE: 0.99765
Corr. Coef.: 0.99945

PIXEL-WISE CALCULATION

Saving and thresholding options

VIEWING PARAMETRIC IMAGE FILES

